

Optical Multiplexer for Laser Phase Locking


Ho Cheong Wei

Supervised by Assoc Prof David Wilkowski and Dr Das Swarup


Introduction

The basic idea of this project is to measure certain wavelengths with a high level of precision. This is done using a self-referenced frequency comb. It references its own frequency to self-generated laser pulses. It doubles the lower frequencies of the comb, so it creates a 'comb' of frequencies with the same frequency difference. Overlapping the comb with a continuous wave laser allows us to phase lock to laser beam to a frequency of the comb. To a comb line, so to speak.

3D Optix Model



Solidworks Model



3DOptix and Solidworks model side by side. 3DOptix was used for laser positions and Solidworks was used for proper mirror positions as well as output.

Details

The setup uses a diffraction grating and mirrors that are meant to be as energy efficient with the wavelengths used as possible. Precautions also have to be made with the diffraction grating so that it does not have unwanted orders of refraction from the grating. In this case, the incoming beam is at an angle of 15 degrees.

Frequencies used -

623nm/679nm/688nm/707nm

Potential Use

High precision measurements are useful for many things. An example would be to use it in a clock.